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PRINCIPIO

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Clark, G. (1981) "Staining Procedures", 4th ed., Williams & Wilkins.

Krafts Woronoff-Dashkoff, Kristine, Clinics in Laboratory Medicine. Vol 22 (2002).

Biological Stains (1977) 9th ed. Ed. by R.D. Lillie; Williams & Wilkins.

CLSI Guidelines and Standards, CLSI, Wayne, P.A

Maree, L.; du Plessis, S.S.; Menkvels, R. and van der Horst, G. Human Reproduction, 25 (6), 1369 – 1382 (2010).

Sousa A. P. M. et al.(2008). Human Reproduction, vol. 24, no. 1

Jimenez-Heffernan J. A. et al., Acta Cytologica, pp. 1-12, (2021)

QUICK PANOPTIC STAIN

PRINCIPLE

The stains which make up the Quick Panoptic stain combine polychromy and the quality of classic haematology staining methods (May-Grünwald, Giemsa, Wright) with a very quick execution time of just 15 seconds. The technique is performed by immersion in the staining solutions.

As with the other Romanowsky stains, the basic dyes bind to the acidic components of cells, nucleic acids, granules in neutrophils and acidic proteins which are stained a relatively deep red-purple colour, while the acid dyes bind to haemoglobin, basic components of cell structures and eosinophil granules. Allowing the differential staining of the different cell types in hematology and cytology.

Furthermore, Quick Panoptic stain allows the polychromatic differentiation of all structures found on sperm cells (acrosome, post-acrosome, vacuoles, midpiece, tail and excess of residual cytoplasm) for a perfect morphological assessment. It can also be used for the differential staining of blood cells and immature germ-cells in sperm samples.

INTENDED USE

Family of in vitro diagnostic medical devices intended for the morphological differentiation of the different cell structures present in hematological samples such as peripheral blood or bone marrow as well as in cytological samples such as biopsies or sperm.

They are qualitative stains that provide information on the physiological or pathological state of the tissue sample when the staining obtained is interpreted by a trained professional. To obtain a diagnosis, the result is utilized alongside other information such as the patient's medical history, physical status as well as other clinical and laboratory data obtained.

For trained and specialized personnel

REAGENTS AND COMPOSITION

Kit 3 x 100 mL

Contains	Ref. 99 30 60
A. Quick panoptic nr. 1	1 x 100 mL
B. Quick panoptic nr. 2	1 x 100 mL
C. Quick panoptic nr. 3	1 x 100 mL

Kit 3 x 500 mL

Contains	Ref. 99 30 90
A. Quick panoptic nr. 1	1 x 500 mL
B. Quick panoptic nr. 2	1 x 500 mL
C. Quick panoptic nr. 3	1 x 500 mL

Quick panoptic Nº1

Composition:	
Fast Green, Cl nº 42 053	<3 mg/L
Methanol	>99 %

Quick panoptic Nº2

Composition:	
Eosine, Cl N. 45 380	<2 g/L
Ethanol 96°	2%
Phosphate Buffer, pH 6,8	0,1M

Quick panoptic Nº3

Composition:	
Methylene Blue Dye, Cl N. 52 015	<1.5 g/L
Azuro B, Cl N. 52 010	<1.5 g/L
Phosphate Buffer, pH 6,8	0,1M

STORAGE AND STABILITY

Reagents which are stored at 15-30°C and protected from light are stable until the expiry date stated on the label. Containers must always be kept tightly closed.

Over time and temperatures below 15°C could lead to the formation of light precipitate in some reagents. This does not affect their functionality. The stain should be filtered before use. It is recommended to bring the reagent to room temperature (15-30°C), agitate and filter before using.

In Use stability should be determined by each user discretion.

MATERIAL NOT SUPPLIED

General-purpose laboratory material.
Slides for blood smears.
Manual or automated staining system.
Microscope.

Stains can be used with a manual staining method or with automatic staining equipment (see QCA stainers available on the website). It is recommended to follow the instructions for use of the automatic equipment supplied by the manufacturer.

PRECAUTIONS

All reagents must be handled with care by trained technicians. Safety instructions are on the product label. The safety data sheet should be consulted before using the reagents. Waste disposal must be carried out according to the local regulations in force.

Any serious incident that has occurred in relation to the device shall be reported to QCA and the competent authority of the member state in which the user is established.

SAMPLE

Air-dried blood and bone marrow smears. Body fluids or tissues smears.

Handling of the samples must be carried out in accordance with the established protocols in each laboratory for the preparation of samples using state-of-the-art technology. It is recommended that they are thin and homogeneous in order to obtain a better fixation of the dye without overstaining.

Staining should be carried out during the two hours following preparation in order to obtain good results. Old smears may stain irregularly.

The ideal sample is capillary blood, but if venous blood is used then EDTA should be used as an anticoagulant. The use of Heparin is not advised.

Handle the samples with care due to their potentially infectious nature.

QUALITY CONTROL

It is recommended to follow the quality control practices defined by the international organisations and periodic staining with routine quality control slide containing tissue processed in a similar manner to the test specimens should be performed.

PREPARATION OF THE WORKING REAGENT

The stain is ready for use, for manual staining and for use in the QCA STAINER.

PROCEDURE

HAEMATOLOGY STAIN

1. Leave to air dry.
2. Set out stains No. 1, No. 2 and No. 3 in Wertheim (or similar) cuvettes.
3. Immerse the rack with the slides for 5 seconds (5 immersions lasting 1 second each) in staining solution No. 1 and drain.
4. Immerse for another 5 seconds (5 immersions lasting 1 second each) in staining solution No. 2 and drain again.
5. Immerse for another 5 seconds (5 immersions lasting 1 second each) in staining solution No. 3.
6. Rinse with tap water and leave to dry.
7. Examine under the microscope.

SPERM CELLS' STAIN

1. Prepare a smear with 15 µl of fresh semen on a standard slide and leave to air dry for at least 10 minutes.
2. Set out stains No. 1, No. 2 and No. 3 in Wertheim (or similar) cuvettes.
3. Immerse the rack with the slides for 5 seconds (5 immersions lasting 1 second each) in staining solution No. 1 and drain.
4. Immerse for another 5 seconds (5 immersions lasting 1 second each) in staining solution No. 2 and drain again.
5. Immerse for another 5 seconds (5 immersions lasting 1 second each) in staining solution No. 3.
6. Rinse with tap water and leave to dry.
7. Examine under the microscope.

3, PG II UN: 1230

H225,H301+H311+H331,H370

P210,P260,P280,P308+P311,P321,P370+P378,P403+P233

GB - QUICK PANOPTIC Nr. 1

Danger

Hazard: Highly flammable liquid and vapour. Toxic if swallowed, in contact with skin or if inhaled. Causes damage to organs.

Precautionary: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Do not breathe dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. If exposed or concerned: Call a POISON CENTER/doctor/... Specific treatment. In case of fire: Use the means described in point 5 of the Safety Data Sheet. Store in a well-ventilated place. Keep container tightly closed.

Contains: methanol

—
GB - QUICK PANOPTIC Nr. 2

—
GB - QUICK PANOPTIC, Nr. 3

For the QCA STAINER, the different methods are found in the system software.

Each user must validate this procedure in his laboratory in order to adjust it to his own standard method, applying the various variants of this procedure, both manual and automatic, allowing him to optimize the results. It must be taken into account that the intensity of the staining is proportional to the staining time.

RESULTS

Red blood cells: Pink or pink-orange.

Platelets: Pale violet or purple.

Neutrophils: Dark violet nucleus. Pink cytoplasm with red-violet granulations.

Eosinophils: Violet nucleus. Blue cytoplasm with red or red-orange granules.

Basophils: Dark blue or purple nucleus. Purple, almost black, granules.

Lymphocytes: Purple nucleus. Sky blue cytoplasm.

Monocytes: Very pale violet nucleus. Sky blue cytoplasm.

Head of sperm cell: Dark violet.

Acrosome of sperm cell: Violet, but a lighter shade of violet than the head.

Mid-piece and tail of the sperm cell: Dark violet.

Background: Pale pink.

NOTES

The result of the staining should be treated as a guide and must be confirmed with additional tests.

The colour results indicated should be treated as a guideline as the intensity and shade of the colour may be influenced by several factors, such as the fixation, the staining time (the staining intensity is proportional to the staining time). The staining intensity can be varied by modifying the number of immersions in stains No. 2 and No. 3, depending on whether it is wished to emphasise the pink or blue shades, the pH value of the staining solution and the buffer solution (if the pH is too basic, the staining will be more blue; and if the pH is too acidic, then the staining will be more pink).

The cuvettes containing the stains, especially that of stain **No. 1**, should be kept covered, in order to prevent excessive evaporation, which could lead to colour deviations compared to normal staining.

Before the content of the cuvette is depleted, add new quantities of the staining solution on a daily basis in order to maintain an adequate level. Renew all the content from time to time. Using buffered water instead of running tap water is recommended for the rinsing processes. Running tap water has a varied pH and salt content, which, along with the high concentrations of chlorine, may alter the results of the staining and produce inconsistent results.

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