

Certificate

Traceability of Manufacturers to the IFCC Reference Measurement Procedure for HbA1c

This certifies that **DiaSys GmbH** using **Innovastar**, participates in the Monitoring Programme to demonstrate traceability. In the Monitoring Programme of 2014 the following performance was seen:

Deviation from IFCC-target	at 30 mmol HbA1c/mol Hb :	1.8
	at 60 mmol HbA1c/mol Hb :	1.0
	at 90 mmol HbA1c/mol Hb :	0.3
Reproducibility, coefficient of variation		1.74%
Linearity, correlation coefficient		0.9974

Date of issue: 10 December 2014

Certification expires: 31 December 2015



IFCC Network Coordinator