

Prämeta tourniquet green



REF: 70.902.00.001 and 70.902.00.002

Product group: Infusion/Transfusion Accessories

Intended use:
peripheral tourniquet

Product description:

The Prämeta tourniquet congests the blood flow of peripheral venous vessels, which is intended to dilate the congested blood vessel and allow better puncture.

REF	Colour	Pack. unit	Individual packing	GTIN	Outer case	GTIN
70.902.00.001	green	10 PCS/pack	10 PCS	04251697615527	10 PCS	04251697615527
70.902.00.002	green	252 PCS/pack	252 PCS	04251697615534	252 PCS	04251697615534

Material

Product	Aluminium alloy, Cotton, POM, Stainless steel, Alu
Package	Cardboard, Corrugated board
Latex	No
DEHP/BBP/DBP	Not applicable

Further information

Regulatory requirements	Not applicable
Notified body	Not applicable
Medical class	I
UMDNS	16-632
eClass	34221214
Supplier	Dahlhausen Medizintechnik GmbH
Disposable product	No
Sterilization	Not applicable
Shelf life	Not applicable
Instruction for use	Yes